
# TB / TUA / TUA VDD / TUS / HELIOS™

Pacing for Cardiology

+++

Temporäre Stimulationskatheter  
*Temporary Transvenous Pacing Catheters*

 **OSYPKA**  
Technology for an active life

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## Temporary Stimulation Catheters

OSYPKA offers a comprehensive array of products for all temporary pacing needs in the fields of interventional adult and pediatric cardiology. Pacing leads are available with varying stiffness. Ideal placement and steerability is achieved through a stylet-driven design in some models. Depending on the catheter type, they are especially suited for:


- Pre- and postoperative stimulation
- Support during pacemaker implantation/exchange
- Rapid pacing during TAVI, valvuloplasty
- Long-term stimulation (max. 30 days) of patients, e.g. with PM pocket infection
- Mapping during electrophysiological procedures
- Reference electrode for 3D Mapping
- Emergency stimulation

OSYPKA temporary transvenous pacing catheters can be used via an OSYPKA extension cable compatible or directly with an OSYPKA pacemakers, such as **PACE 101 H** (single chamber), **PACE 203 H** (dual chamber) and **PACE 300** (triple chamber for biventricular pacing/CRT). A contamination shield (**Cath-Guard®**) can be used with the introducer set (**Adelante® Sigma Plus**), to reduce the risk of contamination.

## TUA II + TUA II VDD Stimulation Catheter

**TUA-II** electrodes are particularly suitable for both recording intracardiac signals and for temporary pacing. The catheters are equipped with a stiff proximal and a flexible distal end for better insertion and positioning. Using a stylet the catheter can easily be positioned at any location inside the heart.

**TUA-II VDD Single Pass** - Quadripolar temporary electrode with fixed connector. With a single catheter, bipolar recording of atrial and ventricular signals as well as bipolar stimulation of the ventricle (VDD mode) is possible.

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### TUS Screw-Tip Stimulation Catheter

**TUS** - Temporary Screw-in Lead with stylet. The distal screw is intended for active fixation in the ventricle as well as in the atrium. TUS is particularly suitable for pacemaker-dependent patients, as well as for patients requiring extended temporary pacing, e.g. in case of pacemaker pocket infections.

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
### TB + TB J Stimulation Catheter

**TB-IIX / TB-IIX-J** - Temporary "All-Rounder" lead with a fixed connector. Straight or with J-Curve for atrial septal wall placement.

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
### HELIOS™ Floating Balloon Stimulation Catheter

Especially convenient in emergencies such as bradyarrhythmia or saltatory conduction diseases. HELIOS™ can be placed in the heart easily and without X-ray control.

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HELIOS™ Convenience Kit also available, please inquire with our product management.

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
As an external stimulation and defibrillation electrode OSYPKA also offers the **R2** electrode for transthoracic pacing with the **OSYPKA PACE 500** and **PACE 500 D** generators.

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