

# R Series® Wi-Fi

## Technical Specifications:



|  |   |
|--|---|
| <b>Operational Bands</b>                               | IEEE 802.11a/b/g/n<br>Dual Band: 2.4 & 5 GHz  |
| <b>Encryption</b>                                      | WPA/TKIP<br>WPA/AES (CCMP)<br>WPA2/TKIP<br>WPA2/AES (CCMP)  |
| <b>Authentication</b>                                  | PSK (Pre-shared Key or "Personal")<br>Enterprise Authentication methods (IEEE 802.1X/RADIUS) as follows: <ul style="list-style-type: none"><li>• EAP-TLS</li><li>• PEAPv0/v1-MSCHAPv2</li></ul> |
| <b>SSID</b>  | Broadcasted by AP<br>Not broadcast (stealth mode)   |
| <b>Preconfigured Profiles<br/>IP Configuration</b>     | Two<br>IP v4<br>Dynamic (external DHCP server)<br>Static IP   |
| <b>Root Certificate Key Length<br/>(if applicable)</b> | 2048-bit is supported – default ZOLL Root Cert uses 2048-bit  |